

H1N1 Puerto Rico Recombinant Protein



RPPB5671

Product Information Protein Information

Product SKU:

RPPB5671

Host:

Baculovirus Insect
Cells.

Protein description:

Recombinant Full-Length H1N1 Puerto Rico 08/1934 is glycosylated with N-linked sugars, produced using baculovirus vectors in insect cells and its Mw is approximately 62 kDa.

Appearance:

Sterile Filtered colorless solution.

Formulation:

The Recombinant H1N1 A/Puerto Rico 08/1934 solution contains 10mM Sodium phosphate, pH 7.2, 150mM NaCl and 0.005% Tween-20.

Purity:

Greater than 90.0% as determined by SDS-PAGE.

Stability:

H1N1 A/Puerto Rico 08/1934 recombinant should be stored at 4°C. Do not freeze!