

# Influenza-B Jilin Recombinant Protein



RPPB5674

---

## Product Information    Protein Information

**Product SKU:**

RPPB5674

**Host:**

Baculovirus Insect  
Cells.

**Protein description:**

Recombinant Full-Length B/Jilin/20/2003 is glycosylated with N-linked sugars, produced using baculovirus vectors in insect cells.

**Appearance:**

Sterile Filtered colorless solution.

**Formulation:**

The Recombinant B/Jilin/20/2003 solution contains 10mM Sodium phosphate, pH 7.4 and 150mM Sodium Chloride.

**Purity:**

Greater than 90.0% as determined by SDS-PAGE.

**Stability:**

B/Jilin/20/2003 Recombinant should be stored at 4°C.