

Influenza-B Influenza B Malaysia Recombinant Protein



RPPB5675

Product Information Protein Information

Product SKU:

RPPB5675

Host:

Baculovirus Insect Cells.

Protein description:

Recombinant Full-Length B/Malaysia/2506/2004 is glycosylated with N-linked sugars, produced using baculovirus vectors in insect cells. The MW is approximately 72,000 Dalton.

Appearance:

Sterile Filtered colorless solution.

Formulation:

The Recombinant B/Malaysia/2506/2004 solution contains 10mM Sodium Phosphate, pH 7.0 and 150mM NaCl, 0.005% Tween-20.

Purity:

Greater than 90.0% as determined by SDS-PAGE.

Stability:

B/Malaysia/2506/2004 Recombinant should be stored at 4°C.