

WNV Pre-M Recombinant Protein



RPPB5677

Product Information Protein Information

Product SKU:

RPPB5677

Host:

Escherichia Coli.

Protein description:

The E.Coli derived 20kda recombinant protein contains the West-Nile N-Terminal Pre-M Virus immunodominant regions. The protein is fused with 6xHis tag at c-terminal.

Formulation:

20mM phosphate buffer pH 7.5.

Purity:

Protein is >95% pure as determined by SDS-PAGE.

Stability:

WNV Pre-M although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.

Amino Acid Sequence:

MVTLSNFQGKVMMTVNATDVTDVITPTAAGKNLCIVRA

MDVGYLCEDTITYECPVLAAGNDPEDIDCWCTKSSVYVRYGRCTKTRHSRRSRRSLTVQTHGESTLANKKGAWLDS
TKATRYLVKTESWILRNPGYALE.