

# Dengue Envelope-1 22kDa Recombinant Protein



RPPB5679

---

## Product Information    Protein Information

**Product SKU:**

RPPB5679

**Protein description:**

The E.coli derived recombinant 22 kDa protein is genetically engineered peptide which is derived from Dengue Type-1 N-terminus Envelope immunodeterminant regions and fused with a 6 aa His Tag. This region also contains a common antigen for Dengue type 1, 2 and 3.

**Formulation:**

Phosphate buffered saline, pH-7.4.

**Purity:**

Protein is >95% pure as determined by 10% PAGE (coomassie staining).

**Stability:**

Dengue Envelope ST1 although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.