

Dengue Envelope-2 22kDa Recombinant Protein



RPPB5680

Product Information Protein Information

Product SKU:

RPPB5680

Host:

Escherichia Coli.

Protein description:

The E.coli derived recombinant 22 kDa protein is genetically engineered peptide which is derived from Dengue Type-2 Envelope. This region also contains a common antigen for Dengue type 2, 3 and 4. This protein is fused to 6xHis tag.

Formulation:

Phosphorous buffered saline, pH-7.4.

Purity:

Protein is >95% pure as determined by 12% PAGE (coomassie staining).

Stability:

Dengue Envelope ST2 although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.