

# Dengue Envelope-2 32kDa Recombinant Protein



RPPB5686

---

## Product Information    Protein Information

**Product SKU:**

RPPB5686

**Host:**

E.coli.

**Protein description:**

Recombinant Dengue Virus Subtype 2 Envelope 32kDa, produced in E.coli, amino acids 45-297 and fused to a 6xHis Tag, containing domains I+II of the dengue envelope.

**Formulation:**

Phosphate buffer.

**Purity:**

Protein is >95% pure as determined by 12% PAGE (coomassie staining).

**Stability:**

Dengue Envelope-2 32kDa although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.