

Dengue Envelope-4 32kDa Recombinant Protein



RPPB5688

Product Information Protein Information

Product SKU:

RPPB5688

Host:

E.coli.

Protein description:

Recombinant Dengue Virus Subtype 4 Envelope 32kDa, produced in E.coli and fused with a 6xHis Tag, containing domains I+II of the dengue envelope.

Formulation:

Phosphate buffered saline, pH-7.4.

Purity:

Protein is >95% pure as determined by 10% PAGE (coomassie staining).

Stability:

Dengue Envelope-4 32kDa although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.