

# Dengue Envelope-1 45kDa Recombinant Protein



RPPB5689

---

## Product Information    Protein Information

**Product SKU:**

RPPB5689

**Host:**

E.coli.

**Protein description:**

Recombinant Dengue Virus Subtype 1 Envelope 45kDa (43-413a.a.) produced in E.coli, contains important epitopes for dengue antibody recognition. This protein is fused to a 6xHis tag.

**Formulation:**

Phosphate buffered saline with 25mM arginine.

**Purity:**

Protein is >95% pure as determined by 10% PAGE (coomassie staining).

**Stability:**

Dengue Envelope-1 45kDa although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.