

Dengue Envelope-2 45kDa Recombinant Protein



RPPB5690

Product Information Protein Information

Product SKU:

RPPB5690

Host:

E.coli.

Protein description:

Recombinant Dengue Virus Subtype 2 Envelope (a.a 43-413) migrates at 45kDa, produced in E.coli, this protein is fused to 6xHis tag at the C-terminus.

Appearance:

Sterile filtered colorless solution.

Formulation:

The protein solution contains Phosphate buffer and 50mM arginine.

Purity:

Protein is >95% pure as determined by 12% PAGE (coomassie staining).

Stability:

Dengue Envelope-2 45kDa although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.