

Dengue Envelope-3 45kDa Recombinant Protein



RPPB5691

Product Information Protein Information

Product SKU:

RPPB5691

Host:

E.coli.

Protein description:

Recombinant Dengue Virus Subtype 3 Envelope 45kDa, produced in E.coli. This protein is fused to 6xHis tag at C-terminus.

Formulation:

Phosphate buffered saline.

Purity:

Protein is >95% pure as determined by 12% PAGE (coomassie staining).

Stability:

Dengue Envelope-3 45kDa although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.