

Dengue Envelope-4 45kDa Recombinant Protein



RPPB5692

Product Information Protein Information

Product SKU:

RPPB5692

Host:

E.coli.

Protein description:

Recombinant Dengue Virus Subtype 4 Envelope 45kDa (a.a 43-413) is produced in E.coli. This protein is fused to a 6xHis tag at C-terminus.

Formulation:

PBS and 25mM arginine.

Purity:

Protein is >95% pure as determined by 12% PAGE (coomassie staining).

Stability:

Dengue Envelope-4 45kDa although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.