

# Dengue Envelope-2 Recombinant Protein



RPPB5699

---

## Product Information    Protein Information

**Product SKU:**

RPPB5699

**Accession:**

ABE02263

**Host:**

Insect cells.

**Protein description:**

Recombinant Dengue Virus Subtype 2 produced in Insect Cells is a polypeptide chain containing amino acids 281-675 and having a molecular weight of approximately 50kDa. Dengue Envelope-2 is purified by proprietary chromatographic technique.

**Appearance:**

Sterile filtered colorless solution.

**Formulation:**

Dengue Envelope-2 protein solution in 1xD-PBS, pH7.4, 0.1% Thimerosal, 5mM EDTA, 1µg/ml of Leupeptin, Aprotinin and Pepstatin A.

**Purity:**

Protein is >95% pure as determined by 12.5% SDS-PAGE.

**Stability:**

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.