

Dengue NS1 c Recombinant Protein



RPPB5704

Product Information Protein Information

Product SKU:

RPPB5704

Host:

Escherichia Coli.

Protein description:

The E.Coli derived recombinant protein contains the NS1 Dengue Virus c-end Type-2 immunodominant regions.

Formulation:

20mM Tris-HCl pH 8.0, 8M urea & 60mM NaCl.

Purity:

Dengue Protein is >90% pure as determined by 10% PAGE (coomassie staining).

Stability:

Dengue although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.

Amino Acid Sequence:

LWSNGVLESE MVIKKNFAGP KSQHNNRPGY HTQTAGPWHL GKLEMDDFDC EGTTVVVTEG CGNRGPSLRT
TTASGKLITE WCCRSCTLPP LRYRGEDGCW YGMEIRPLKE KEENLVSSLV TAHHHHHH