

# Dengue NS1 n Recombinant Protein



RPPB5705

---

## Product Information    Protein Information

**Product SKU:**

RPPB5705

**Host:**

Escherichia Coli.

**Protein description:**

The E.coli derived recombinant protein contains the NS1 Dengue Virus n-end Type-2 immunodominant regions. The Recombinant Dengue Virus NS1n is fused to a GST tag.

**Formulation:**

100mM Tris-Hcl pH 8.0, 60mM NaCl, and 5% Glycerol.

**Purity:**

Protein is >90% pure as determined by SDS-PAGE.

**Stability:**

Dengue although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.