

Dengue NS1 ST1 Recombinant Protein



RPPB5707

Product Information Protein Information

Product SKU:

RPPB5707

Accession:

ACA48761

Host:

Insect cells.

Protein description:

Recombinant Dengue Virus NS1 Subtype 1 produced in Insect Cells is a polypeptide chain containing amino acids 777-1131 and having a molecular weight of approximately 50kDa. Dengue NS1 ST1 is purified by proprietary chromatographic technique.

Appearance:

Sterile filtered colorless solution.

Formulation:

Dengue NS1 ST1 protein solution (1mg/ml) in 1xPBS, pH7.4, 0.1% Thimerosal, 5mM EDTA, 1µg/ml of Leupeptin, Aprotinin and Pepstatin A.

Purity:

Protein is >95% pure as determined by 12.5% SDS-PAGE.

Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.