

Dengue NS1 Recombinant Protein



RPPB5709

Product Information Protein Information

Product SKU:

RPPB5709

Host:

Escherichia Coli.

Protein description:

The E.Coli derived recombinant protein contains the NS1 Dengue Virus full length Type-3 immunodominant regions, migrates as 45kDa on 12% SDS-PAGE gel. The dengue protein is fused to a 6xHis tag at C-terminus.

Formulation:

Phosphate buffered saline and 50mM arginine.

Purity:

Protein is >90% pure as determined by 12% PAGE (coomassie staining).

Stability:

Dengue NS1, ST3 although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.