

# Dengue C-terminal Recombinant Protein



RPPB5713

---

## Product Information    Protein Information

**Product SKU:**

RPPB5713

**Host:**

E.coli.

**Protein description:**

The Recombinant Dengue 2 Envelope Hydrophobic Domain e.coli derived contains 225 amino acids (228-485) from dengue 2 envelope. The protein is fused to a His tag at C-terminal and purified by standard chromatography techniques.

**Appearance:**

Sterile Filtered clear solution.

**Formulation:**

PBS.

**Purity:**

Protein is >95% pure as determined by 10% PAGE (coomassie staining).

**Stability:**

Dengue 2 Envelope Recombinant although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.