

Human Procalcitonin Recombinant Protein



RPPB5748

Product Information Protein Information

Product SKU:

RPPB5748

Host:

Escherichia Coli.

Protein description:

Procalcitonin Human Recombinant produced in E.Coli is a single, non-glycosylated, polypeptide chain containing 114 amino acids fragment (3-116) having a molecular mass of 17.13 kDa and an amino-terminal hexahistidine tag. The PCT is purified by standard chromatographic techniques.

Appearance:

Sterile Filtered clear solution.

Synonyms:

Procalcitonin, PCT.

Formulation:

PCT is supplied in 20mM Tris-HCl pH 8.0 and 50% glycerol.

Purity:

Greater than 95.0% as determined by SDS-PAGE.

Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. Please prevent freeze-thaw cycles.