

Mouse CALB2 Recombinant Protein



RPPB5865

Product Information Protein Information

Product SKU:

RPPB5865

Accession:

Q08331

Host:

Escherichia Coli.

Protein description:

Calretinin Mouse Recombinant full length protein expressed in E.coli, shows a 57 kDa band on SDS-PAGE. The Calretinin is purified by proprietary chromatographic techniques.

Appearance:

Sterile Filtered clear solution.

Synonyms:

Calretinin, CR, Calb2, calbindin 2.

Formulation:

Calretinin protein at 100µg/ml in 50mM Tris-HCl, pH7.5 and 10mM L-glutathione (reduced).

Stability:

Store vial at -20°C to -80°C. When stored at the recommended temperature, this protein is stable for 12 months. Please prevent freeze-thaw cycles.