

# Polcalcin Phl p 7 Recombinant Protein



RPPB5943

---

## Product Information    Protein Information

**Product SKU:**

RPPB5943

**Accession:**

O82040

**Host:**

Sf9 insect cells.

**Protein description:**

Recombinant Polcalcin Phl p 7 produced in SF9 is a glycosylated, polypeptide chain having a calculated molecular mass of 8,805 Dalton. Polcalcin Phl p 7 is purified by proprietary chromatographic techniques.

**Appearance:**

Sterile Filtered clear solution.

**Synonyms:**

Polcalcin Phl p 7, Calcium-binding pollen allergen Phl p 7, P7, Phl p 7.

**Formulation:**

Polcalcin Phl p 7 is supplied in 20mM HEPES buffer pH-7.9 and 6M Urea.

**Purity:**

Greater than 80.0% as determined by SDS-PAGE.

**Stability:**

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. Avoid multiple freeze-thaw cycles.