

HIV-2 gp36 Recombinant Protein



RPPB5947

Product Information Protein Information

Product SKU:

RPPB5947

Host:

Synthetic.

Protein description:

HIV-2 gp36 is full length chemically synthesized polypeptide sequence of HIV-2 envelope immunodominant regions.

Appearance:

Sterile filtered colorless clear solution.

Formulation:

(1mg/1ml) in H₂O.

Purity:

Greater than 95.0% as determined by HPLC.

Stability:

HIV-2 gp36 although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.