

Influenza-B Qingdao Recombinant Protein



RPPB5949

Product Information Protein Information

Product SKU:

RPPB5949

Protein description:

Allantoic fluid of 10 days old embryonated eggs, inoculated with influenza B virus, strain B/Qingdao/102/91. The Influenza B Virus was purified by Ultracentrifugation with 10-40 % sucrose gradient.

Appearance:

Sterile Filtered colorless solution.

Formulation:

The B/Qingdao/102/91 solution contains 0.1M NaCl, 10mM Tris HCl, 1mM EDTA pH-8, 0.1% sodium azide (NaN₃) and 0.005% thimerosal.

Purity:

Greater than 90.0% as determined by Analysis by SDS-PAGE.

Stability:

B/Qingdao/102/91 although stable 4°C for 4 weeks, should be stored below -18°C. Please prevent freeze-thaw cycles.