

Influenza-B Malaysia Recombinant Protein



RPPB5952

Product Information Protein Information

Product SKU:

RPPB5952

Protein description:

Allantoic fluid of 10 days old embryonated eggs, inoculated with influenza B virus, strain B/Malaysia/2506/04. The Influenza B Virus was purified by Ultracentrifugation with 10-40 % sucrose gradient.

Appearance:

Sterile Filtered colorless solution.

Formulation:

The B/Malaysia/2506/04 solution contains STE, 0.1% sodium azide (NaN₃) and 0.005% thimerosal.

Purity:

Greater than 90.0% as determined by SDS-PAGE.

Stability:

B/Malaysia/2506/04 although stable 4°C for 4 weeks, should be stored below -18°C. Please prevent freeze-thaw cycles.