

H1N1 Beijing Recombinant Protein



RPPB5955

Product Information Protein Information

Product SKU:

RPPB5955

Protein description:

Allantoic fluid of 10 days old embryonated eggs, inoculated with influenza A virus, strain A/Beijing/262/95. The Influenza Virus was purified by Ultracentrifugation with 10-40 % sucrose gradient.

Appearance:

Opaque suspension.

Formulation:

The H1N1 A/Beijing/262/95 solution contains STE, 0.09% sodium azide (NaN₃) and 0.005% thimerosal. Stir thoroughly prior to its use.

Purity:

Greater than 90.0% as determined by SDS-PAGE.

Stability:

H1N1 A/Beijing/262/95 although stable 4°C for 4 weeks, should be stored below -18°C. Please prevent freeze-thaw cycles.