

H1N1 Taiwan Recombinant Protein



RPPB5957

Product Information Protein Information

Product SKU:

RPPB5957

Protein description:

Allantoic fluid of 10 days old embryonated eggs, inoculated with influenza A virus, strain A/Taiwan/1/86. The Influenza Virus was purified by Ultracentrifugation with 10-40 % sucrose gradient.

Appearance:

Sterile Filtered colorless solution.

Formulation:

The H1N1A/Taiwan/1/86 solution contains STE, 0.1% sodium azide (NaN₃) and 0.005% thimerosal.

Purity:

Greater than 90.0% as determined by SDS-PAGE.

Stability:

A/Taiwan/1/86 although stable 4°C for 4 weeks, should be stored below -18°C. Please prevent freeze-thaw cycles.