

Influenza-A Solomon Islands Recombinant Protein



RPPB5958

Product Information Protein Information

Product SKU:

RPPB5958

Protein description:

Allantoic fluid of 10 days old embryonated eggs, inoculated with influenza A virus, strain A/Solomon Islands/03/06. The Influenza Virus was purified by Ultracentrifugation with 10-40 % sucrose gradient.

Appearance:

Sterile Filtered colorless solution.

Formulation:

The H1N1 A/Solomon Islands/03/06 solution contains 0.1M NaCl, 10mM Tris-HCl, 1mM EDTA pH-8 (STE), 0.1% sodium azide (NaN₃) and 0.005% thimerosal.

Purity:

Greater than 90.0% as determined by SDS-PAGE.

Stability:

A/Solomon Islands/03/06 although stable 4°C for 4 weeks, should be stored below -18°C. Please prevent freeze-thaw cycles.