

H3N2 Shandong Recombinant Protein



RPPB5959

Product Information Protein Information

Product SKU:

RPPB5959

Protein description:

Allantoic fluid of 10 days old embryonated eggs, inoculated with influenza A virus, strain A/Shandong/9/93. The Influenza Virus was purified by Ultracentrifugation with 10-40 % sucrose gradient.

Appearance:

Opaque suspension.

Formulation:

The H3N2A/ Shandong/9/93 solution (1.0mg/ml) contains STE, 0.1% sodium azide (NaN₃) and 0.005% thimerosal. Stir thoroughly prior to its use.

Purity:

Greater than 90.0%.

Stability:

A/Shandong/9/93 although stable 4°C for 4 weeks, should be stored below -18°C. Please prevent freeze-thaw cycles.