

H3N2 Kiev Recombinant Protein



RPPB5960

Product Information Protein Information

Product SKU:

RPPB5960

Protein description:

Allantoic fluid of 10 days old embryonated eggs, inoculated with influenza A virus, strain A/Kiev/301/94 like /Johannesburg/33/94. The Influenza Virus was purified by Ultracentrifugation with 10-40 % sucrose gradient.

Appearance:

Sterile Filtered colorless solution.

Formulation:

The H3N2 A/Kiev/301/94 solution contains 0.1M NaCl, 10mM Tris-Hcl, 1mM EDTA pH-8, 0.1% sodium azide (NaN₃) and 0.005% thimerosal.

Purity:

Greater than 90.0% as determined by Analysis by SDS-PAGE.

Stability:

A/Kiev/301/94 although stable 4°C for 4 weeks, should be stored below -18°C. Please prevent freeze-thaw cycles.