

# H3N2 Panama Recombinant Protein



RPPB5961

---

## Product Information    Protein Information

**Product SKU:**

RPPB5961

**Protein description:**

Allantoic fluid of 10 days old embryonated eggs, inoculated with influenza A virus, strain A/Panama/2007/99. The Influenza Virus was purified by Ultracentrifugation with 10-40 % sucrose gradient.

**Appearance:**

Opaque suspension.

**Formulation:**

The H3N2 A/Panama/2007/99 solution contains STE, 0.09% sodium azide (NaN<sub>3</sub>) and 0.005% thimerosal. Stir thoroughly prior to its use.

**Purity:**

Greater than 90.0% as determined by Analysis by SDS-PAGE.

**Stability:**

A/Panama/2007/99 although stable 4°C for 4 weeks, should be stored below -18°C. Please prevent freeze-thaw cycles.