

# Parainfluenza Type-2 Recombinant Protein



RPPB5965

---

## Product Information    Protein Information

**Product SKU:**

RPPB5965

**Protein description:**

MA 104 cells inoculated with parainfluenza 2 virus, strain II ALTB cc 2056.

**Appearance:**

Sterile Filtered colorless solution.

**Formulation:**

The Parainfluenza type-2 solution contains STE, 0.02% sodium azide (NaN<sub>3</sub>).

**Purity:**

Greater than 90.0% as determined by Analysis by SDS-PAGE.

**Stability:**

Parainfluenza type-2 although stable 4°C for 4 weeks, should be stored below -18°C. Please prevent freeze-thaw cycles.