

# Parainfluenza Type-3 Recombinant Protein



RPPB5966

---

## Product Information   Protein Information

**Product SKU:**

RPPB5966

**Protein description:**

MA 104 cells inoculated with Parainfluenza-3 virus, strain III v2932.

**Appearance:**

Opaque (non-clear) solution. Do not centrifuge! Mix thoroughly before use!

**Formulation:**

The Parainfluenza type-3 solution contains STE, 0.09% sodium azide (NaN<sub>3</sub>) and 0.005% thimerosal. Do not centrifuge! Mix thoroughly before use!

**Purity:**

Greater than 90.0% as determined by Analysis by SDS-PAGE.

**Stability:**

Parainfluenza type-3 although stable 4°C for 4 weeks, Should be stored below -18°C. Please prevent freeze-thaw cycles.