

HDPT0651

---

## Product Information

### Uniprot:

P19338

### Expression Host:

HEK293

### Target

NCL

### Usage:

Research use only

### Background:

Nucleolin (NCL), a eukaryotic nucleolar phosphoprotein, is involved in the synthesis and maturation of ribosomes. It is located mainly in dense fibrillar regions of the nucleolus. Human NCL gene consists of 14 exons with 13 introns and spans approximately 11kb. The intron 11 of the NCL gene encodes a small nucleolar RNA, termed U20.

### Description:

Recombinant human NCL protein with N-terminal 6xHis tag

### Synonyms:

C23;Nsr1

### Tag:

N-6xHis tag

### Molecular Weight:

The protein has a predicted molecular mass of 77.3 kDa after removal of the signal peptide. The apparent molecular mass of His-NCL is approximately 100-130 kDa due to glycosylation.

### Purity:

The purity of the protein is greater than 85% as determined by SDS-PAGE and Coomassie blue staining.

### Molecular Characterization:

6xHis tag NCL(Val2-Glu710)

### Formulation:


Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization.

### Storage:

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.

## Product Images

---



**Figure 1.** Human NCL Protein, His Tag on SDS-PAGE under reducing condition.