

Purified Anti-Human CD35 Antibody [E11]

AGEL0257

Description

This Purified Anti-Human CD35 Antibody [E11] is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGEL0257
Contents:	100µg, 25µg Bradford Reagent: 1 vial (2ml)
Category:	Monoclonal Antibody
Clonality:	Monoclonal
Clone:	E11
Synonyms:	C3BR, C3b/C4b receptor, CD35, CR1, Complement receptor type 1
Applications:	FCM
Reactivity:	Human
Immunogen:	Recombinant Human CD35 protein

Antibody Data

Uniprot ID:	-
Gene ID:	-
Swissprot:	P17927
Host Species:	Mouse
Isotype:	Mouse IgG1, κ
Isotype Control:	-

Conjugation:	-
Conjugation Information:	-
Buffer:	Phosphate-buffered solution, pH 7.2, containing 0.05% non-protein stabilizer. Dialyze to completely remove the stabilizer prior to labeling.
Purification:	>98%, Protein A/G purified
Target:	CD35
Cellular Localization:	Cell membrane
Tissue Specificity:	-
Verified Samples:	Verified Samples in FCM: Human peripheral blood granulocytes cell
Concentration:	≥ 1 mg/mL

Preparation & Storage

Storage: Store at 4°C valid for 12 months or -20°C valid for long term storage, avoid freeze / thaw cycles.

Shipping: Ice bag

Recommended Dilution: FCM 2 µg/mL(0.5×10⁶-1×10⁶ cells)

Recommended Usage:

Application	Recommended Usage
-	-

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Notes: Centrifuge before opening to ensure complete recovery of vial contents.