

Anti-Mouse CD120b Antibody [TR75-54.7] -Low Endotoxin Purified

AGEL2165

Description

This Anti-Mouse CD120b Antibody [TR75-54.7] -Low Endotoxin Purified is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGEL2165
Contents:	100µg, 25µg Bradford Reagent: 1 vial (2ml)
Category:	Monoclonal Antibody
Clonality:	Monoclonal
Clone:	TR75-54.7
Synonyms:	CD120b, TNF-R2, TNF-RII, TNFR-II, Tnfrsf1b, Tumor necrosis factor receptor 2, Tumor necrosis factor receptor superfamily member 1B, p75, p80 TNF-alpha receptor
Applications:	FCM Block
Reactivity:	Mouse
Immunogen:	-

Antibody Data

Uniprot ID:	-
Gene ID:	-
Swissprot:	P25119
Host Species:	Armenian Hamster
Isotype:	Armenian Hamster IgG

Isotype Control:	-
Conjugation:	-
Conjugation Information:	-
Buffer:	Sterile PBS, pH 7.2. < 1.0 EU per mg of the antibody as determined by the LAL method .
Purification:	-
Target:	CD120b/TNFR2
Cellular Localization:	-
Tissue Specificity:	-
Verified Samples:	-
Concentration:	≥ 1 mg/mL

Preparation & Storage

Storage: Store at 4°C valid for 12 months or -20°C valid for long term storage, avoid freeze / thaw cycles. This preparation contains no preservatives, thus it should be handled under aseptic conditions.

Shipping: Ice bag

Recommended Dilution: -

Recommended Usage:

Application	Recommended Usage
-	-

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Notes: Centrifuge before opening to ensure complete recovery of vial contents.