

Purified Armenian Hamster IgG Isotype Control [PIP]

AGEL2698

Description

This Purified Armenian Hamster IgG Isotype Control [PIP] is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

| | |
|----------------------|---|
| SKU: | AGEL2698 |
| Contents: | 100µg, 25µg Bradford Reagent: 1 vial (2ml) |
| Category: | Monoclonal Antibody |
| Clonality: | Monoclonal |
| Clone: | PIP |
| Synonyms: | - |
| Applications: | FCM |
| Reactivity: | - |
| Immunogen: | - |

Antibody Data

| | |
|-------------------------|----------------------|
| Uniprot ID: | - |
| Gene ID: | - |
| Swissprot: | - |
| Host Species: | Armenian Hamster |
| Isotype: | Armenian Hamster IgG |
| Isotype Control: | - |

| | |
|---------------------------------|--|
| Conjugation: | - |
| Conjugation Information: | - |
| Buffer: | Phosphate-buffered solution, pH 7.2, containing 0.05% non-protein stabilizer. Dialyze to completely remove the stabilizer prior to labeling. |
| Purification: | - |
| Target: | - |
| Cellular Localization: | - |
| Tissue Specificity: | - |
| Verified Samples: | - |
| Concentration: | ≥ 1 mg/mL |

Preparation & Storage

Storage: Store at 4°C valid for 12 months or -20°C valid for long term storage, avoid freeze / thaw cycles.

Shipping: Ice bag

Recommended Dilution: -

| | | |
|---------------------------|--------------------|--------------------------|
| Recommended Usage: | Application | Recommended Usage |
| | - | - |

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Notes: Centrifuge before opening to ensure complete recovery of vial contents.