

Anti-Human/Mouse KLRG-1 Antibody [2F1] -Low Endotoxin, Azide Free

AGEL4999

Description

This Anti-Human/Mouse KLRG-1 Antibody [2F1] -Low Endotoxin, Azide Free is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGEL4999
Contents:	25µg, 100µg, 1mg Bradford Reagent: 1 vial (2ml)
Category:	Monoclonal Antibody
Clonality:	Monoclonal
Clone:	2F1
Synonyms:	Klrg, CLEC15A, MAFA, MAFAL, KLRG1, HNK-1, Leu-7, Klrg, Mafa, Klrg1
Applications:	FCM
Reactivity:	Human,Mouse
Immunogen:	Recombinant Mouse KLRG-1 protein

Antibody Data

Uniprot ID:	-
Gene ID:	-
Swissprot:	Q96E93O88713
Host Species:	Syrian Hamster
Isotype:	Syrian Hamster IgG
Isotype Control:	-

Conjugation:	-
Conjugation Information:	-
Buffer:	Sterile PBS, pH 7.2. < 1.0 EU per mg of the antibody as determined by the LAL method.
Purification:	>98%, Protein A/G purified
Target:	KLRG-1
Cellular Localization:	-
Tissue Specificity:	-
Verified Samples:	-
Concentration:	≥ 1 mg/mL

Preparation & Storage

Storage: Store at 4°C valid for 12 months or -20°C valid for long term storage, avoid freeze / thaw cycles. This preparation contains no preservatives, thus it should be handled under aseptic conditions.

Shipping: Ice bag

Recommended Dilution: FCM 2 µg/mL(0.5×10⁶-1×10⁶ cells)

Recommended Usage:

Application	Recommended Usage
-	-

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Notes: Centrifuge before opening to ensure complete recovery of vial contents.