

Anti-Human CD154 Antibody [24-31] -Low Endotoxin, Azide Free

AGEL5029

Description

This Anti-Human CD154 Antibody [24-31] -Low Endotoxin, Azide Free is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGEL5029
Contents:	25µg, 100µg, 1mg Bradford Reagent: 1 vial (2ml)
Category:	Monoclonal Antibody
Clonality:	Monoclonal
Clone:	24-31
Synonyms:	IMD, HIGM, T-Cell Antigen Gp, Tumor Necrosis Factor Ligand Superfamily Member, TNFSF, CD154, CD40L, HIGM1, IGM, IMD3, T-BAM, TNFSF5, TRAP, gp39, hCD40L, CD40LG, CD40 Ligand, CD40-L, T-Cell Antigen Gp39, TNF-Related Activation Protein, Tumor Necrosis Factor Ligand Superfamily Member 5, T-BAM, gp39, TRAP, TNFSF5, TNFSF5, CD40LG
Applications:	FCM
Reactivity:	Human
Immunogen:	Recombinant Human CD154 protein

Antibody Data

Uniprot ID:	-
Gene ID:	-
Swissprot:	P29965
Host Species:	Mouse

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Isotype:	Mouse IgG1, κ
Isotype Control:	-
Conjugation:	-
Conjugation Information:	-
Buffer:	Sterile PBS, pH 7.2. < 1.0 EU per mg of the antibody as determined by the LAL method .
Purification:	>98%, Protein A/G purified
Target:	CD154
Cellular Localization:	-
Tissue Specificity:	-
Verified Samples:	Verified Samples in FCM: Human peripheral blood lymphocytes
Concentration:	≥ 1 mg/mL

Preparation & Storage

Storage: Store at 4°C valid for 12 months or -20°C valid for long term storage, avoid freeze / thaw cycles. This preparation contains no preservatives, thus it should be handled under aseptic conditions.

Shipping: Ice bag

Recommended Dilution: FCM 2 µg/mL(0.5×10⁶-1×10⁶ cells)

Recommended Usage:

Application	Recommended Usage
-	-

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Notes: Centrifuge before opening to ensure complete recovery of vial contents.