

Purified Anti-Human CD1a Antibody [Ab571]

AGEL5071

Description

This Purified Anti-Human CD1a Antibody [Ab571] is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGEL5071
Contents:	25µg, 100µg, 1mg Bradford Reagent: 1 vial (2ml)
Category:	Monoclonal Antibody
Clonality:	Monoclonal
Clone:	Ab571
Synonyms:	FCB, HTA, T-cell surface antigen T6/Leu, CD1A, CD1, FCB6, HTA1, R4, T6, CD 1a, T-cell surface glycoprotein CD1a, T-cell surface antigen T6/Leu-6, T6, HTA1, a polypeptide, CD1A Antigen, CD1a molecule, cluster of differentiation 1 A, cortical thymocyte antigen CD1A, differentiation antigen CD1 alpha 3, epidermal dendritic cell marker CD1a, FCB 6, HTA 1, hTa1 thymocyte antigen, OTTHUMP00000018907, R 4, T 6
Applications:	FCM
Reactivity:	Human
Immunogen:	Recombinant Human CD1a protein

Antibody Data

Uniprot ID:	-
Gene ID:	-
Swissprot:	P06126

Host Species:	Mouse
Isotype:	Mouse IgG1, κ
Isotype Control:	-
Conjugation:	-
Conjugation Information:	-
Buffer:	Phosphate-buffered solution, pH 7.2, containing 0.05% non-protein stabilizer. Dialyze to completely remove the stabilizer prior to labeling.
Purification:	>98%, Protein A/G purified
Target:	CD1a
Cellular Localization:	-
Tissue Specificity:	-
Verified Samples:	Verified Samples in FCM: Molt-4
Concentration:	≥ 1 mg/mL

Preparation & Storage

Storage: Store at 4°C valid for 12 months or -20°C valid for long term storage, avoid freeze / thaw cycles.

Shipping: Ice bag

Recommended Dilution: FCM 2 µg/mL(0.5×10⁶-1×10⁶ cells)

Recommended Usage:

Application	Recommended Usage
-	-

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Notes: Centrifuge before opening to ensure complete recovery of vial contents.