

# Purified Anti-Human CD95 Antibody [DX2]

AGEL5098

## Description

---

This Purified Anti-Human CD95 Antibody [DX2] is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

## Product Information

---

<b>SKU:</b>	AGEL5098
<b>Contents:</b>	25µg, 100µg, 1mg Bradford Reagent: 1 vial (2ml)
<b>Category:</b>	Monoclonal Antibody
<b>Clonality:</b>	Monoclonal
<b>Clone:</b>	DX2
<b>Synonyms:</b>	APT, TNF, CD antigen CD, Tumor necrosis factor receptor superfamily member, TNFRSF, APO, ALPS1A, APO-1, APT1, CD95, FAS1, FASTM, TNFRSF6, FAS, APO1, FASR, TNFRSF6 and FAS, TNFR6, Tumor necrosis factor receptor superfamily member 6, Apo-1 antigen, Apoptosis-mediating surface antigen FAS, FASLG receptor, CD antigen CD95, APO-1, APT1, FASTM, TNFRSF6, ALPS 1A, APO 1, Apo 1 antigen, APO 1 cell surface antigen, Apo1 antigen, APO1 cell surface antigen, Apoptosis antigen 1, Apoptosis mediating surface antigen FAS, APT 1, CD 95, CD 95 antigen, CD95 antigen, Delta Fas, Delta Fas/APO 1/CD95, Delta Fas/APO1/CD95, Fas (TNF receptor superfamily, FAS 1, FAS 827dupA, Fas AMA, FAS Antigen, Fas cell surface death receptor, FAS/CD95, member 6, member 6), sFAS, Surface antigen APO1, TNF receptor superfamily, TNFRSF 6
<b>Applications:</b>	<b>FCM</b>
<b>Reactivity:</b>	Human
<b>Immunogen:</b>	Recombinant Human CD95 protein

## Antibody Data

---

<b>Uniprot ID:</b>	-
<b>Gene ID:</b>	-
<b>Swissprot:</b>	P25445
<b>Host Species:</b>	Mouse
<b>Isotype:</b>	Mouse IgG1, κ
<b>Isotype Control:</b>	-
<b>Conjugation:</b>	-
<b>Conjugation Information:</b>	-
<b>Buffer:</b>	Phosphate-buffered solution, pH 7.2, containing 0.05% non-protein stabilizer. Dialyze to completely remove the stabilizer prior to labeling.
<b>Purification:</b>	>98%, Protein A/G purified
<b>Target:</b>	CD95
<b>Cellular Localization:</b>	-
<b>Tissue Specificity:</b>	-
<b>Verified Samples:</b>	Verified Samples in FCM: Human peripheral blood lymphocytes
<b>Concentration:</b>	≥ 1 mg/mL

## Preparation & Storage

---

<b>Storage:</b>	Store at 4°C valid for 12 months or -20°C valid for long term storage, avoid freeze / thaw cycles.
<b>Shipping:</b>	Ice bag
<b>Recommended Dilution:</b>	FCM 2 µg/mL(0.5×10 <sup>6</sup> -1×10 <sup>6</sup> cells)

**Recommended Usage:**

Application	Recommended Usage
-	-

**Protein  
Quantification  
(Optional):**

To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

**Notes:**

Centrifuge before opening to ensure complete recovery of vial contents.