

Purified Anti-Mouse CD11c Antibody [N418]

AGEL5146

Description

This Purified Anti-Mouse CD11c Antibody [N418] is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGEL5146
Contents:	25µg, 100µg, 1mg Bradford Reagent: 1 vial (2ml)
Category:	Monoclonal Antibody
Clonality:	Monoclonal
Clone:	N418
Synonyms:	Leukocyte adhesion receptor p150+, CD11 antigen-like family member C, CD11c, Integrin alpha-X, Leukocyte adhesion receptor p150+95, Itgax, p150-95, CR4, integrin &, alpha-X-X
Applications:	FCM
Reactivity:	Mouse
Immunogen:	Recombinant Mouse CD11c protein

Antibody Data

Uniprot ID:	-
Gene ID:	-
Swissprot:	Q9QXH4
Host Species:	Armenian Hamster
Isotype:	Armenian Hamster IgG

Isotype Control:	-
Conjugation:	-
Conjugation Information:	-
Buffer:	Phosphate-buffered solution, pH 7.2, containing 0.05% non-protein stabilizer. Dialyze to completely remove the stabilizer prior to labeling.
Purification:	>98%, Protein A/G purified
Target:	CD11c
Cellular Localization:	-
Tissue Specificity:	-
Verified Samples:	Verified Samples in FCM: C57/BL6 Mouse splenocytes
Concentration:	≥ 1 mg/mL

Preparation & Storage

Storage: Store at 4°C valid for 12 months or -20°C valid for long term storage, avoid freeze / thaw cycles.

Shipping: Ice bag

Recommended Dilution: FCM 2 µg/mL(0.5×10⁶-1×10⁶ cells)

Recommended Usage:

Application	Recommended Usage
-	-

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Notes: Centrifuge before opening to ensure complete recovery of vial contents.