

Purified Anti-Mouse CD172a Antibody [P84]

AGEL5161

Description

This Purified Anti-Mouse CD172a Antibody [P84] is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGEL5161
Contents:	25µg, 100µg, 1mg Bradford Reagent: 1 vial (2ml)
Category:	Monoclonal Antibody
Clonality:	Monoclonal
Clone:	P84
Synonyms:	Myd, SHPS, Ptpns, Bit, Myd1, Ptpns1, Shps1, Sirp, Sirpa, SIRP&, beta-1-1, SIRP&, alpha, SIRP&, beta, A1835480, CD172a, P84, SHP-1, SHPS-1
Applications:	FCM
Reactivity:	Mouse
Immunogen:	Recombinant Mouse CD172a protein

Antibody Data

Uniprot ID:	-
Gene ID:	-
Swissprot:	P97797
Host Species:	Rat
Isotype:	Rat IgG1, κ

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Isotype Control:	-
Conjugation:	-
Conjugation Information:	-
Buffer:	Phosphate-buffered solution, pH 7.2, containing 0.05% non-protein stabilizer. Dialyze to completely remove the stabilizer prior to labeling.
Purification:	>98%, Protein A/G purified
Target:	CD172a
Cellular Localization:	-
Tissue Specificity:	-
Verified Samples:	Verified Samples in FCM: C57/BL6 Mouse bone marrow
Concentration:	≥ 1 mg/mL

Preparation & Storage

Storage:	Store at 4°C valid for 12 months or -20°C valid for long term storage, avoid freeze / thaw cycles.					
Shipping:	Ice bag					
Recommended Dilution:	FCM 2 µg/mL(0.5×10 ⁶ -1×10 ⁶ cells)					
Recommended Usage:	<table><tr><th>Application</th><th>Recommended Usage</th></tr><tr><td>-</td><td>-</td></tr></table>		Application	Recommended Usage	-	-
Application	Recommended Usage					
-	-					
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol					
Notes:	Centrifuge before opening to ensure complete recovery of vial contents.					