

# Anti-Human CD105 Antibody [SN6] -Low Endotoxin, Azide Free

AGEL5175

## Description

---

This Anti-Human CD105 Antibody [SN6] -Low Endotoxin, Azide Free is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

## Product Information

---

<b>SKU:</b>	AGEL5175
<b>Contents:</b>	25µg, 100µg Bradford Reagent: 1 vial (2ml)
<b>Category:</b>	Monoclonal Antibody
<b>Clonality:</b>	Monoclonal
<b>Clone:</b>	SN6
<b>Synonyms:</b>	HHT, RP11 228B, Osler Rendu Weber syndrome, END, HHT1, ORW1, CD105, ENG, endoglin, AI528660, AI662476, CD 105, CD105 antigen, EGLN, FLJ41744, ORW, Osler Rendu Weber syndrome 1, RP11 228B15.2, S-endoglin, SN6, AI528660, AI662476, CD 105, CD105, CD105 antigen, EGLN, END, Endoglin, Eng, FLJ41744, HHT1, ORW, ORW1, Osler Rendu Weber syndrome 1, RP11 228B15.2, S-endoglin, SN6, ENG, END
<b>Applications:</b>	<a href="#">FCM</a>
<b>Reactivity:</b>	Human
<b>Immunogen:</b>	Recombinant Human CD105 protein

## Antibody Data

---

<b>Uniprot ID:</b>	-
<b>Gene ID:</b>	-
<b>Swissprot:</b>	P17813

<b>Host Species:</b>	Mouse
<b>Isotype:</b>	Mouse IgG1, $\kappa$
<b>Isotype Control:</b>	-
<b>Conjugation:</b>	-
<b>Conjugation Information:</b>	-
<b>Buffer:</b>	Sterile PBS, pH 7.2. < 1.0 EU per mg of the antibody as determined by the LAL method.
<b>Purification:</b>	>98%, Protein A/G purified
<b>Target:</b>	CD105
<b>Cellular Localization:</b>	-
<b>Tissue Specificity:</b>	-
<b>Verified Samples:</b>	Verified Samples in FCM: U937
<b>Concentration:</b>	$\geq 1$ mg/mL

## Preparation & Storage

---

**Storage:** Store at 4°C valid for 12 months or -20°C valid for long term storage, avoid freeze / thaw cycles. This preparation contains no preservatives, thus it should be handled under aseptic conditions.

**Shipping:** Ice bag

**Recommended Dilution:** FCM 2  $\mu$ g/mL( $0.5 \times 10^6$ - $1 \times 10^6$  cells)

**Recommended Usage:**

Application	Recommended Usage
-	-

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

**Notes:** Centrifuge before opening to ensure complete recovery of vial contents.