

Purified Anti-Mouse Ly-49C/I Antibody [5E6]

AGEL5273

Description

This Purified Anti-Mouse Ly-49C/I Antibody [5E6] is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGEL5273
Contents:	25µg, 100µg, 1mg Bradford Reagent: 1 vial (2ml)
Category:	Monoclonal Antibody
Clonality:	Monoclonal
Clone:	5E6
Synonyms:	Klra, Ly-49c, Ly49C, Klra3
Applications:	FCM
Reactivity:	Mouse
Immunogen:	Recombinant Mouse Ly-49C/I protein

Antibody Data

Uniprot ID:	-
Gene ID:	-
Swissprot:	Q64329
Host Species:	Mouse
Isotype:	Mouse IgG2a, κ
Isotype Control:	-

Conjugation:	-
Conjugation Information:	-
Buffer:	Phosphate-buffered solution, pH 7.2, containing 0.05% non-protein stabilizer. Dialyze to completely remove the stabilizer prior to labeling.
Purification:	>98%, Protein A/G purified
Target:	Ly-49C/I
Cellular Localization:	-
Tissue Specificity:	-
Verified Samples:	-
Concentration:	≥ 1 mg/mL

Preparation & Storage

Storage: Store at 4°C valid for 12 months or -20°C valid for long term storage, avoid freeze / thaw cycles.

Shipping: Ice bag

Recommended Dilution: -

Recommended Usage:	Application	Recommended Usage
	-	-

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Notes: Centrifuge before opening to ensure complete recovery of vial contents.