

Purified Anti-Mouse CD105 Antibody [MJ7/18]

AGEL5285

Description

This Purified Anti-Mouse CD105 Antibody [MJ7/18] is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGEL5285
Contents:	25µg, 100µg, 1mg Bradford Reagent: 1 vial (2ml)
Category:	Monoclonal Antibody
Clonality:	Monoclonal
Clone:	MJ7/18
Synonyms:	HHT, ORW1, AI528660, AI662476, CD105, END, Endo, Endoglin, HHT1, ORW, S-endoglin, Edg, Eng, END, HHT1, ORW, ORW1, CD 105, CD105 antigen, EGLN, FLJ41744, Osler Rendu Weber syndrome 1, RP11 228B15.2, SN6
Applications:	FCM
Reactivity:	Mouse
Immunogen:	Recombinant Human CD105 protein

Antibody Data

Uniprot ID:	-
Gene ID:	-
Swissprot:	Q63961
Host Species:	Rat
Isotype:	Rat IgG2a, κ

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Isotype Control:	-
Conjugation:	-
Conjugation Information:	-
Buffer:	Phosphate-buffered solution, pH 7.2, containing 0.05% non-protein stabilizer. Dialyze to completely remove the stabilizer prior to labeling.
Purification:	>98%, Protein A/G purified
Target:	CD105
Cellular Localization:	-
Tissue Specificity:	-
Verified Samples:	Verified Samples in FCM: Bend.3
Concentration:	≥ 1 mg/mL

Preparation & Storage

Storage:	Store at 4°C valid for 12 months or -20°C valid for long term storage, avoid freeze / thaw cycles.				
Shipping:	Ice bag				
Recommended Dilution:	FCM 2 µg/mL(0.5×10 ⁶ -1×10 ⁶ cells)				
Recommended Usage:	<table border="1"> <thead> <tr> <th>Application</th><th>Recommended Usage</th></tr> </thead> <tbody> <tr> <td>-</td><td>-</td></tr> </tbody> </table>	Application	Recommended Usage	-	-
Application	Recommended Usage				
-	-				
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol				
Notes:	Centrifuge before opening to ensure complete recovery of vial contents.				