

Anti-CHMP2B [R06-4U6] Monoclonal Antibody

AGMB00077

Description

This Anti-CHMP2B [R06-4U6] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB00077
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	ALS17, CHMP2b, DMT1, hVps2-2, VPS2B, CGI-84, CHMP2.5, chromatin modifying protein 2B, Chromatin-modifying protein 2b
Applications:	WB ICC/IF FC IP
Research Area:	Signal Transduction
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat
Clone:	R06-4U6
Clonality:	Monoclonal Antibody
SwissProt ID:	Q9UQN3
Immunogen:	A synthesized peptide derived from human CHMP2B
Gene ID:	25978

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: CHMP2B

Host Species: Rabbit

Isotype: IgG

Purification: Affinity Purified

Conjugated: Unconjugated

Modification: Unmodified

Molecular Weight: Calculated MW: 24 kDa, Observed MW: 28 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 10mM PBS, pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000
	IF	1:50-1:200
	IP	1:50
	FC	1:50-1:100

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.