

Anti-NXF1 [R05-1C3] Monoclonal Antibody

AGMB00136

Description

This Anti-NXF1 [R05-1C3] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB00136
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	DmNXF1, MEX67, Mex67p, Mvb1, NXF1, Sbr, TAP, DKFZp667O0311, DmNXF1, MEX67,
Applications:	WB IHC-P ICC/IF FC IP
Research Area:	Tags & Cell Markers
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	R05-1C3
Clonality:	Monoclonal Antibody
SwissProt ID:	Q9UBU9
Immunogen:	A synthetic peptide of human NXF1
Gene ID:	10482

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: NXF1
Host Species: Rabbit
Isotype: IgG
Purification: Affinity Purified
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: Calculated MW: 70 kDa, Observed MW: 70 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:1000-1:5000
	ICC/IF	1:200-1:500
	FC	1:20-1:50
	IP	1:10-1:100
	IHC	1:50-1:100

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.