

Anti-DBF4 [R06-6I2] Monoclonal Antibody

AGMB00137

Description

This Anti-DBF4 [R06-6I2] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB00137
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	ASK, DBF4A, ZDBF1, Protein DBF4 homolog A, Activator of S phase kinase, Chiffon homolog A, DBF4-type zinc finger-containing protein 1
Applications:	WB
Research Area:	Cell Biology
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	R06-6I2
Clonality:	Monoclonal Antibody
SwissProt ID:	Q9UBU7
Immunogen:	A synthetic peptide of human DBF4
Gene ID:	10926

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: DBF4
Host Species: Rabbit
Isotype: IgG
Purification: Affinity Purified
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: Calculated MW: 77 kDa, Observed MW: 77 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:2000-1:5000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.