

# Anti-CGBP [R06-4K1] Monoclonal Antibody

## AGMB00164

### Description

---

This Anti-CGBP [R06-4K1] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Product Information

---

|                       |                                                                                                        |
|-----------------------|--------------------------------------------------------------------------------------------------------|
| <b>SKU:</b>           | AGMB00164                                                                                              |
| <b>Contents:</b>      | 20µl, 50µl, 100µl<br>Bradford Reagent: 1 vial (2ml)                                                    |
| <b>Synonyms:</b>      | CFP1, CGBP, CpG binding protein, CXXC finger 1 (PHD domain), CXXC1, hCGBP, HsT2645, PCCX1, PHF18, SPP1 |
| <b>Applications:</b>  | <b>WB</b> <b>IHC-F</b> <b>IHC-P</b> <b>ICC/IF</b> <b>FC</b> <b>IP</b>                                  |
| <b>Research Area:</b> | Epigenetics                                                                                            |
| <b>Form:</b>          | Liquid                                                                                                 |

### Antibody Data

---

|                      |                             |
|----------------------|-----------------------------|
| <b>Reactivity:</b>   | Human, Mouse, Rat           |
| <b>Clone:</b>        | R06-4K1                     |
| <b>Clonality:</b>    | Monoclonal Antibody         |
| <b>SwissProt ID:</b> | Q9P0U4                      |
| <b>Immunogen:</b>    | A synthetic peptide of CGBP |
| <b>Gene ID:</b>      | 30827                       |

#### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

|                          |                                            |
|--------------------------|--------------------------------------------|
| <b>Gene Name:</b>        | CXXC1                                      |
| <b>Host Species:</b>     | Rabbit                                     |
| <b>Isotype:</b>          | IgG                                        |
| <b>Purification:</b>     | Affinity Purified                          |
| <b>Conjugated:</b>       | Unconjugated                               |
| <b>Modification:</b>     | Unmodified                                 |
| <b>Molecular Weight:</b> | Calculated MW: 76 kDa, Observed MW: 76 kDa |

## Preparation & Storage

---

**Storage:** Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.  
Store Bradford Reagent at Room Temperature for 1 Year.

**Storage Buffer:** Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

| <b>Antibody Dilution Ratio:</b>           | Application                                                                                                                                                                                                                                             | Antibody Dilution Ratio |
|-------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------|
|                                           | WB                                                                                                                                                                                                                                                      | 1:500-1:1000            |
|                                           | IHC-P                                                                                                                                                                                                                                                   | 1:200-1:1000            |
|                                           | ICC/IF                                                                                                                                                                                                                                                  | 1:100-1:200             |
|                                           | FC                                                                                                                                                                                                                                                      | 1:50-1:100              |
|                                           | IP                                                                                                                                                                                                                                                      | 1:20-1:50               |
| <b>Protein Quantification (Optional):</b> | To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <a href="https://www.assaygenie.com/bradford-protein-assay-protocol/">https://www.assaygenie.com/bradford-protein-assay-protocol/</a> to view the full protocol. |                         |