

Anti-ACAT2/Acetyl-CoA acetyltransferase [R88-4C-2] Monoclonal Antibody

AGMB00311

Description

This Anti-ACAT2/Acetyl-CoA acetyltransferase [R88-4C-2] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB00311
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	Acat2, ACTL, Acetyl-CoA transferase-like protein, Cytosolic acetoacetyl-CoA thiolase, Acat2, acetoacetyl CoA thiolase, acetoacetyl Coenzyme A thiolase,
Applications:	WB IHC-P
Research Area:	Cardiovascular
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat
Clone:	R88-4C-2
Clonality:	Monoclonal Antibody
SwissProt ID:	Q9BWD1
Immunogen:	A synthesized peptide derived from human ACAT2

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene ID: 39

Gene Name: ACAT2

Host Species: Rabbit

Isotype: IgG

Purification: Affinity Chromatography

Conjugated: Unconjugated

Modification: Unmodified

Molecular Weight: Calculated MW: 41 kDa, Observed MW: 38kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 10mM PBS, pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

Antibody Dilution Ratio:

Application	Antibody Dilution Ratio
WB	1:500-1:1000
IHC	1:100-1:200

Protein Quantification (Optional):

To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.